

SKIN CLEANSING

# Cimolin

Universal wash gel for gentle skin cleaning



**DR.SCHNELL**  
PROFESSIONAL SKIN CARE

## PRODUCT DESCRIPTION

CIMOLIN is a wash gel with mild cleaning power and excellent skin-compatibility. CIMOLIN is economical in use and has a pleasant scent. CIMOLIN cleans skin intensely yet gently with a skin-neutral pH-value and creamy foam. CIMOLIN is free of alkali and soap and contains moisturizing substances.

## AREA OF APPLICATION

CIMOLIN wash gel is ideal for light to middle grade soiling and is suitable for all areas of working life. CIMOLIN is especially skin-compatible and is perfect for damaged, dry skin.

## PRODUCT ADVANTAGES

- Mild and pleasant
- Dermatologically tested
- pH skin neutral
- Soap- and alkali-free
- Fresh scent
- Helps skinregeneration
- Available in a practical 500 ml bottle with flip top cap as individual bottle for 500 ml Euro-wall dispenser or with a hand-pump dispenser
- Ecological 1000 ml dispenser bottle
- Hygienic, easy and germ-free dosing with the new 1000 ml DR.SCHNELL wall dispenser

## PRODUCT INFORMATION

- Appearance: light blue, viscous liquid
- pH-value skin neutral 5.5

## SUGGESTED APPLICATION

Wash dirty skin with approx. 2 ml of CIMOLIN and a very little amount of water, then add more water and continue to wash skin until dirt is removed. Thereafter rinse with a lot of water and dry skin thoroughly. After washing, do not forget to apply moisturizing lotion!

## INGREDIENTS

Aqua, Sodium Laureth Sulfate, Sodium Chloride, Cocoamidopropyl Betaine, Coco-Glucoside, Potassium Lactate, Sodium Laureth Glutamate, Parfum, Benzyl Alcohol (and) 2-Bromo-2-nitropropane-1,3-diol (and) Iodopropynyl Butylcarbamate, Hexyl Cinnamal, Limonene, C.I. Direct Blue 199

## PRODUCT SAFETY AND STORAGE

CIMOLIN has been thoroughly tested on human skin in practical, working use and results show excellent skin-compatibility. CIMOLIN falls under the cosmetic directive and meets the requirements of the LFGB. CIMOLIN can be stored for up to 30 months unopened and at room temperature.

## DELIVERY SIZES AND PACKAGING

500 ml bottle	PU: 10x500 ml	Item No. 50180
1000 ml dispenser bottle*	PU: 6x1000 ml	Item No. 50185
5 litre canister		Item No. 50200

\*Can only be used with the 1000 ml DR.SCHNELL wall dispenser

## ACCESSORIES

DR.SCHNELL wall dispenser for 1000 ml CIMOLIN-bottles	Item No. 80044
Dispenser (hand pump) for 500 ml CIMOLIN-bottles	Item No. 94787

Our Information and recommendations represent today's standard and our knowledge and experience and for which we cannot be held liable. The user/employee must take necessary precaution prior to using the products with regards to application possibilities and area of application.  
Printed: 11/2007.